Sterilizer cycle log

complyday.com

Office: Sterilizer (make, model, serial): Month and year:

Date	Cycle or load no.	Cycle type	Contents	Temp	Time at temp	Pressure	Internal indicator OK (Y/N)	External indicator OK (Y/N)	Initials	Notes

Check the gauges, display or printout every load. Don't use a package if the cycle didn't reach the manufacturer's settings or its indicator didn't change color. Source: CDC, Best Practices for Sterilization Monitoring in Dental Settings.

Weekly spore test log

complyday.com

Office: Sterilizer (make, model, serial): Year:

Date run	Cycle or load no.	Spore test lot no.	Control lot no.	Mail-in or in-office	Date of result	Test result (neg/pos)	Control result	Initials	Action if positive

Run a spore test at least weekly on each sterilizer, and on every load with an implantable device. If a test is positive, take the sterilizer out of service and repeat the test right away on the same cycle. Source: CDC.

Sterilizer maintenance log

complyday.com

Office: Sterilizer (make, model, serial): Year:

Date	Task (from the manufacturer's instructions)	Done by	Result or parts replaced	Next due	Initials

Use the tasks and intervals in your sterilizer's instructions for use, for example chamber cleaning, door gasket checks, water reservoir changes and filter cleaning. Keep service and repair records with this log.